



DB089: Glu-Glu (E6)

Background:

Glu-Glu is an epitope tag composed of a 6 residue peptide, EYMPME, derived from the Polyoma middle T protein.

Origin:

Glu-Glu (E6) is provided as an affinity purified rabbit polyclonal antibody, raised against a peptide derived from the Polyoma middle T protein.

Product Details:

Each vial contains 200 µg/ml of affinity purified rabbit IgG, Glu-Glu (E6) DB089, in 1 ml PBS containing 0.1 % sodium azide and 0.2% gelatin.

Specificity:

Glu-Glu (E6) is recommended for use by western blotting and immunoprecipitation. The suggested western blotting starting dilution is at 1:400.

Storage:

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.