



DB092: T7 (E11)

Background:

T7 is an epitope tag composed of an 11 residue peptide, MASMTGGQMG, derived from the T7 major capsid protein of the T7 phage.

Origin:

T7 (E11) is provided as an affinity purified rabbit polyclonal antibody, raised against a peptide derived from the T7 major capsid protein of the T7 phage.

Product Details:

Each vial contains 200 µg/ml of affinity purified rabbit IgG, T7 (E11) DB092, in 1 ml PBS containing 0.1 % sodium azide and 0.2% gelatin.

Specificity:

T7 (E11) is recommended for use by western blotting and immunoprecipitation. The suggested western blotting starting dilution is at 1:400.

Storage:

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.