



## ***DB093: V5 (E14)***

### **Background:**

V5 is an epitope tag composed of a 14 residue peptide, GKPIPPLLGLDST, derived from the P/V proteins of paramyxovirus SV5.

### **Origin:**

V5 (E14) is provided as an affinity purified rabbit polyclonal antibody, raised against a peptide derived from the P/V proteins of paramyxovirus SV5.

### **Product Details:**

Each vial contains 200 µg/ml of affinity purified rabbit IgG, V5 (E14) DB093, in 1 ml PBS containing 0.1 % sodium azide and 0.2% gelatin.

### **Specificity:**

V5 (E14) is recommended for use by western blotting and immunoprecipitation. The suggested western blotting starting dilution is at 1:1000.

### **Storage:**

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.