



DB121: Protein-A (FL7)

Background:

Protein A is a cell wall constituent of *Staphylococcus aureus*. Protein A is a 42 K protein that has a high affinity for the Fc portion of immunoglobulins, IgG in particular.

Origin:

Protein-A (FL7) is provided as an affinity purified chicken IgG raised against purified recombinant protein A.

Product Details:

Each vial contains 200 µg/ml of affinity purified chicken IgG polyclonal, Protein-A (FL7) DB121, in 1 ml PBS containing 0.1 % sodium azide and 0.2% gelatin.

Specificity:

Protein-A (FL7) DB121 recognizes *Staphylococcus aureus* protein A

Use:

Protein-A (FL7) DB121 is recommended for use by western blotting, immunoprecipitation and ELISA. The suggested western blotting starting dilution is at 1:1000.

Storage:

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.

LABORATORY REAGENT FOR IN-VITRO RESEARCH USE ONLY