



## ***DB2002: anti-Mouse IgG+HRP***

### **Product:**

HRP conjugated affinity purified Goat anti-Mouse IgG h+l.

### **Origin:**

Goats were immunized with highly purified Mouse IgG.  
Antibodies were affinity purified and conjugated with horseradish peroxidase (HRP).

### **Product Details:**

Each vial contains 200 µg in 0.2ml of affinity purified anti-Mouse IgG, in a buffer containing 50mM HEPES pH 7.0, 0.1 M sodium chloride, 0.2% BSA and 0.01 % Thimerosal.

**Use:** Recommended starting dilution for Western blotting 1:5,000-1:20,000.

### **Storage:**

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.