



DB2005: anti-Chicken IgY+HRP

Product:

HRP conjugated affinity purified Goat anti-Chicken IgY (IgG) h+l.

Origin:

Goats were immunized with highly purified Chicken serum IgY (IgG).
Antibodies were affinity purified and conjugated with horseradish peroxidase (HRP).

Product Details:

Each vial contains 200 µg in 0.2ml of affinity purified anti-Chicken IgY (IgG), , in a buffer containing 50mM HEPES pH 7.0, 0.1 M sodium chloride, 0.2% BSA and 0.01 % Thimerosal.

Use:

Recommended starting dilution for Western blotting 1:5,000-1:20,000. Will detect IgY (IgG) in Chicken serum and egg yolks. It also reacts with the papain produced fab and fc fragments plus isolated heavy (IgY only) and light chains.

Storage:

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.