



### ***DB2006: anti-Goat IgG+AP***

**Product:**

AP conjugated affinity purified Rabbit anti-Goat IgG h+l.

**Origin:**

Rabbits were immunized with highly purified Goat IgG.

Antibodies were affinity purified and conjugated with alkaline phosphatase (AP).

**Product Details:**

Each vial contains 200 µg in 0.2ml of affinity purified anti-Goat IgG, in a buffer containing 50mM HEPES pH 7.0, 0.1 M sodium chloride, 0.2% BSA and 0.01 % Thimerosal.

**Use:** Recommended starting dilution for Western blotting 1:5000-1:20,000.

**Storage:**

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.