



DB2009: anti-Rabbit IgG+AP

Product:

AP conjugated affinity purified Goat anti-Rabbit IgG h+l.

Origin:

Goats were immunized with highly purified Rabbit IgG.
Antibodies were affinity purified and conjugated with alkaline phosphatase (AP).

Product Details:

Each vial contains 200 µg in 0.2ml of affinity purified anti-Rabbit IgG, in a buffer containing 50mM HEPES pH 7.0, 0.1 M sodium chloride, 0.2% BSA and 0.01 % Thimerosal.

Use: Recommended starting dilution for Western blotting 1:5,000-1:20,000.

Storage:

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.